



Tubular film with arc and port seal closures

BAG DETAILS

NUMBER OF BAG PORTS: 1 2 3

BAG HEIGHT (mm)

BAG WIDTH (mm)

BAG VOLUME (mL)

BAG VOLUME - MAX (mL)

PRIMARY APPLICATION

- Cryogenic storage
- Biobanking
- Tissue banking
- Cell freezing
- Transport
- Research
- Other

OPERATING TEMP

- over 34 °C
- 34 °C to Ambient
- Ambient to -80 °C
- 80 °C to -200 °C

ADDITIONAL DETAILS

SEND SAMPLE BAGS TO

NAME

COMPANY

EMAIL

PHONE

ADDRESS

PRIMARY BAG ID:

CLP -

PRINT FORM

SAVE FORM